

## Experiment 10

### AIM



To prepare a solution of common salt of 10% composition by mass.

### THEORY



The concentration of a solution is the amount of solute present in a given amount (mass or volume) of the solution. Mass percent concentration is defined as a mass (g) of the solute per 100 g mass of the solution. A 10% solution by mass means, 10 g of solute dissolved in 90 g of solvent to result into 100 g of solution.

### MATERIALS REQUIRED



Common salt, distilled water, watch glass, stirring rod, physical balance, measuring cylinder (100 mL), and a beaker (250 mL).

### PROCEDURE



1. Calculate the volume of solvent (distilled water) and mass of solute (sodium chloride or common salt) required to prepare 100 g of 10% by mass solution. This may be done as follows–

10 g solute is required for 100 g solution (distilled water + salt). Thus the amount of water required would be  $100\text{ g} - 10\text{ g} = 90\text{ g}$ . Since the density of distilled water is  $1\text{ g/mL}$ , therefore the volume of distilled water (required to prepare 100 g of 10% by mass solution with 10 g of common salt) is 90 mL.

2. Weigh an empty watch glass on a physical balance. Also weigh 10 g of sodium chloride (common salt) on the watch glass now.
3. Take 90 mL (90 g) of distilled water in a 250 mL beaker with the help of a measuring cylinder.
4. Transfer 10 g of salt from the watch glass to the beaker containing 90 mL distilled water.
5. Stir the solution until all the salt dissolve in it.

## OBSERVATIONS AND CALCULATIONS



- |  |           |
|--|-----------|
| (i) Mass of empty watch glass ( $m_1$ )                    | = _____ g |
| (ii) Mass of watch glass + sodium chloride ( $m_1 + 10$ g) | = _____ g |
| (iii) Mass of sodium chloride (common salt)                | = 10 g.   |

## RESULTS AND DISCUSSION



The concentration of prepared solution is 10% by mass of common salt in water. This way of expressing concentration is one amongst many.

## PRECAUTIONS AND SOURCES OF ERROR



- Use of physical balance must be done with all precautions. Ask your teacher acquaint you with the working of a physical balance.
- The readings of the physical balance should be taken only when its pointer comes to rest.
- The measuring cylinder should be placed on a horizontal surface while measuring the volume of the distilled water and solution.

## NOTE FOR THE TEACHER

- In place of common salt, some students may be suggested to perform this experiment with sugar.
- This experiments requires to use a physical balance to weigh 10 g of common salt or sugar using a watch glass. A physical balance is a sophisticated equipment. It is suggested to acquaint students with the working of a physical balance. They may be asked to practice with this balance before performing this experiment.
- In case, if distilled water is not available, the experiment may be performed with filtered water or drinking water. Its density may be assumed as 1 g/mL at the experimental temperature.

## QUESTIONS

- Why should the density of a 10% common salt be more than the density of pure water at a specified temperature? Offer qualitative explanation.
- A student is asked to prepare 250 mL sugar solution 15% by mass concentration. How much amount of sugar and water should be taken for preparation of the solution?
- If 50 mL of water is added to the above solution. What will be the change in the mass percentage of the solute?
- 830 g of salt solution contains 50 g of common salt in it. Calculate its concentration in terms of mass percentage?