EXPERIMENT 2

A_{IM}

To determine the resistance of a given wire using a metre bridge and hence determine the resistivity of the material of the wire.

APPARATUS AND MATERIAL REQUIRED

Metre bridge, a wire about 1m long (of material whose specific resistance is to be determined), a resistance box, a rheostat, galvanometer, a jockey, one-way key, a cell or battery eliminator, thick connecting wires, sand paper, screw gauge.

DESCRIPTION OF APPARATUS

Metre bridge

It consists of one metre long constantan wire AC of uniform cross-sectional area mounted on a wooden board with a scale (Fig. E 2.1). The two ends of the wire are attached to terminals A and C. Thick metal strips bent at right angles are used to provide two gaps E and F to connect resistors forming a Wheatstone's




Fig. E 2.1 A metre bridge

bridge (Fig. E 2.2). The terminal B between the gaps is used for connecting galvanometer and other end of the galvanometer is connected to a jockey J.

PRINCIPLE

A metre bridge works on the principle of Wheatstone's bridge. As shown in Fig. E 2.2, it consists of four resistors P,Q,R and S connected in the form of a network ABCD. The terminals A and C are connected to two terminals of a cell through a key K_1 . Terminals B and D are connected to a sensitive galvanometer G through a key K_2 .

If there is no deflection in the galvanometer G, then balance condition for Wheatstone's bridge is

$$\frac{P}{Q} = \frac{R}{S}$$
 (E 2.1)

We use this relation (E 2.1) to determine S if P, Q and R are known.

The unknown resistance S is connected in the gap E and a resistance box (R_{Box}) in gap F of the metre bridge. The terminal B is connected to one terminal of the galvanometer G. The other terminal of the galvanometer is connected to a jockey J which slides along the wire AC. A source of dc current is connected between A and C through a key K_1 so as to provide a constant potential drop along AC (Fig. E 2.1).

A resistor (or wire) of known resistance is inserted in the gap F by taking out corresponding key from the resistance box R_{Box} . The jockey is moved on the wire AC to obtain a condition of no-deflection in the galvanometer. It happens when the jockey is kept at a point D called the null point. In this condition;

Fig.E 2.2 The Wheatstone's bridge

$$\frac{P}{Q} = \frac{R}{S} = \frac{Resistance\ of\ wire\ of\ length\ DC}{Resistance\ of\ wire\ of\ length\ AD}$$

(E 2.2)

Unknown resistance *S* of the wire, having uniform cross-sectional area, is then given by

$$S = R \times \frac{l}{100 - l}$$
 (E 2.3)

Because, for a wire of uniform cross-sectional area, resistance is proportional to length.

Thus, knowing l and R, and using Eq. (E 2.3), the unknown resistance S can be determined.

Resistivity

The specific resistance or resistivity ho of the material of the given wire is

$$\rho = \frac{Sa}{L}$$

where *S* is the resistance of the wire of length *L* and $\alpha = \pi r^2$, (*r* being the radius) is the area of cross-section.

PROCEDURE

- 1. Find the average diameter of the wire with a screw gauge. From this, obtain the value of its radius r.
- 2. Clean the insulation at the ends of connecting wires with a piece of sand paper. Tighten all plugs of the resistance box (R_{BOX}) by pressing each plug.
- 3. Set up the circuit as shown in Fig. E 2.1 with unknown resistance wire of known length in gap E.
- 4. Next, introduce some resistance *R* in the circuit from the resistance box. Bring the jockey J in contact with terminal A first and then with terminal C. Note the direction in which pointer of the galvanometer gets deflected in each case. Make sure that jockey remains in contact with the wire for a fraction of a second. If the galvanometer shows deflection on both sides of its zero mark for these two points of contact of the jockey, null point will be somewhere on the wire AC. If it is not so, adjust resistance *R* so that the null point is somewhere in the middle of the wire AC, say, between 30 cm and 70 cm.
- 5. If there is one-sided deflection, check the circuit again, especially junctions, for their continuity.
- 6. Repeat step 4 for four different values of resistance *R*.
- 7. Interchange the position of the resistances *S* and *R* and repeat steps 4 to 6 for the same five values of *R*. While interchanging *S* and *R*, ensure that the same length of wire of resistance *S* is now in the gap F. The interchange takes care of unaccounted resistance offered by terminals.

OBSERVATIONS

- 1. Length of the wire of unknown resistance, L = ...cm
- 2. Measurement of diameter of wire of unknown resistance

Least count of the screw gauge (L.C.) = ... mm

Zero error of the screw gauge = ...mm

Zero correction of the screw gauge = ...mm

Table E 2.1: Diameter of Wire

S	Re	eading along direction	one	Read perpe	Mean Diameter		
	Main scale reading <i>p</i> (mm)	Coinciding circular scale division n	Diameter $d_1 = p + n \times L.C.$ (mm)	scale	Coinciding circular scale division n'	Diameter $d_2 = p' + n' \times L.C.$ (mm)	$d = \frac{d_1 + d_2}{2}$ (mm)
1 2 3							

Mean diameter (corrected for zero error) = ... mm Radius of wire r = ... mm

Table E 2.2: Unknown resistance

Sl. No.	S in the left gap					S in the right gap									
	(Ho) Resistance R	Position of Balance point D	Balancing length AD (l)	$ \begin{array}{ccc} \widehat{\mathbb{B}} & \text{Length} \\ \widehat{\mathbb{B}} & \text{DC} = 100 - l \end{array} $	$(0.5) S_1 = R \times \frac{l}{100 - l}$	Position of balance point D'	Balancing Belancth AD' =l'	$ \stackrel{\bigcirc}{\mathbb{B}} \text{D'C'} = 100 - l' $	$\frac{\partial}{\partial B} S_2 = R \times \frac{100 - l'}{l'}$	Mean resistance $S = \frac{S_1 + S_2}{2}$	ρ Ω m	ΔS_1 Ω	$\Delta S_{_{2}}$ Ω	$\Delta S = \Omega$	$\Delta \rho \Omega \mathrm{m} \Omega$
1			•												
2															
5															
	-								Mean						

31

CALCULATIONS

$$L = \dots \text{ cm}$$
 $r = \dots \text{ mm}$ $S = \dots \Omega$

Substituting these values in $\rho = S \frac{\pi r^2}{L}$, calculate the value of ρ in ohm-m.

Error

$$\frac{\Delta \rho}{\rho} = \frac{\Delta S}{S} + \frac{2\Delta r}{r} + \frac{\Delta L}{L}$$

Errors Δr , ΔL are the least counts of the measuring instruments and error ΔS is the maximum of the values obtained by following equations:

$$\Delta S_1 = \left[\frac{\Delta l}{l} + \frac{\Delta l}{(100 - l)} \right] S_1$$

$$\Delta S_2 = \left[\frac{\Delta l'}{l'} + \frac{\Delta l'}{(100 - l')}\right] S_2$$

If it is presumed that resistance offered by the resistance box is the same as that indicated on it, then $\Delta R = 0$

Therefore, maximum error is $\Delta S = \Delta S_1 + \Delta S_2$

RESULT

- 1. The unknown resistance of the given wire is found to be $S + \Delta S = ... \pm ... \Omega$
- 2. The resistivity of the material of the wire is $\rho \pm \Delta \rho = \dots \Omega m$ Here S and ρ are mean values. ΔS and $\Delta \rho$ are maximum of the five values of error.

RECAUTIONS

- 1. All the connections and plugs should be tight.
- 2. Jockey should be moved gently over the metre bridge wire.
- 3. The plug in the key (K_1) should be inserted only at the time of taking observations.
- 4. Null points should be in the middle of the wire (30 cm to 70 cm).

Sources of Error

- 1. The metre bridge wire may not be of uniform area of cross-section.
- 2. Effect of end resistances due to copper strips, connecting screws, may affect the measurement.
- 3. The length L of the wire should not include the lengths below the terminals when placed in gap E or F.
- 4. The resistances of end pieces/metal strips may not be negligible. The error introduced by it can be reduced by interchanging the known and unknown resistances in gaps E and F.
- 5. The length measurements *l* and *l'* may have error if the metre bridge wire is not taut and along the scale in the metre bridge.
- 6. Galvanometer pointer is expected to be at zero when no current flows through it. However, many times it is observed that it is not so. In such cases, pointer has to be adjusted to zero by gently moving the screw below the scale with the help of a screw driver. Otherwise null point must be obtained by tapping the jockey on the wire.

DISCUSSION

- 1. It is generally preferred that *R* and *S* combination should be so chosen that the balance point comes near the middle of the metre bridge wire. Why? Is the sensitivity of detecting balance point best when *R* and *S* have same order of resistance?
- 2. The accuracy with which the null point can be detected also depends upon the sensitivity of the galvanometer used. To investigate this, find the distance through which the jockey has to be moved to cause a just perceptible deflection of the galvanometer. How can the error due to this be minimised? Does it have something to do with the potential difference between A and C?
- 3. The length l, that is, the position of the jockey, may be read upto, \pm 0.1 cm (say). How much uncertainty will be introduced in the result due to this error?
- 4. It is advised that key should be unplugged when readings are not being taken to avoid unnecessary heating of the wire. Why? In what way, will the heating affect the null point? Will it be significant?

SELF ASSESSMENT

- 1. The bridge wire may not be exactly 100 cm long, so record its exact length and use this in the calculation.
- 2. If the metre bridge wire is not of uniform area of cross section, how will it affect the observations?




Fig. E 2.3 Defective scale of a metre bridge

- 3. If the metre bridge wire is not made of a material of uniform density, how will it affect the observations?
- 4. If the same experiment is performed with AC = 50 cm wire instead of 1 m, what changes would be there in the result?
- 5. The scale attached along the wire may not indicate the exact length as shown in Fig. E 2.3. How would you minimise the error?
- 6. Why is the metre bridge not suitable for measuring very high/very low resistances?
- 7. Why is the metre bridge preferred for determining the resistance of a wire over the Ohm's law circuit?

SUGGESTED ADDITIONAL EXPERIMENTS/ACTIVITIES

- 1. Determine the resistance of an electric bulb of given voltage and power. Attach crocodile clips at the end of two copper wires. Use these to connect the bulb in the gap E as shown in Fig. E 2.4. Why does the resistance obtained by you differ from the theoretical value?
- 2. Making use of your observations plot a graph between (100 l)/l and R. Find the slope of the graph. What does it represent?
- 3. The resistance of a wire changes with temperature. You may study this by using a metre bridge. You may use a beaker with mustard oil, a suitable thermometer (0-300°C) and resistance wire of known length. The experimental set-up shown in Fig. E 2.5 may be used.




Fig. E 2.4




Fig. E 2.5